

Dräger DrugTest[®] 5000

Oral Fluid Analysis System for Detecting Drugs

- Most sensitive oral fluid POCT device for THC: 5 ng/ml
- Clearly measures results without interpretation
- Easy and simple sample collection
- Integrated data management
- Proven Dräger Quality
- Training and support from SBL's team of experts



Fast, Accurate, and Controlled Oral Fluid Drug Screening

The DrugTest 5000, manufactured by Dräger Safety, is a mobile state-of-the-art opto-electronic oral fluid drug screening device. It is fast, reliable, and entirely eliminates any possibility of human error in reading or interpreting results.

The DrugTest 5000 is the most sensitive oral fluid POCT device available globally and is the first approved device in the Criminal Code of Canada by the Justice Department for roadside drug screening.

Benefits

- Screens for recent vs. historical use
- No known adulterates or sample substitutions
- Witnessed collection with results in minutes
- Less invasive with no dignity or gender issues
- Free online training and certification
- Integral laboratory confirmation and MRO process
- Collector turns blue when sufficient sample is obtained



Data Management



Each test result can be printed multiple times for donor and administrator. The device stores 500 test results for PC download.


The Benefits of Detecting Recent Use

It is important to recognize the difference associated with testing for drugs-of-abuse using oral fluid vs. urine. Oral fluid, as a specimen type, can detect drug use within *minutes of consumption*, whereas urine requires the ingested substances to be metabolized by the liver before detection. This is a process that can take hours.

The ability to detect recent use means employers can determine fit-for-duty much more quickly and easily. The detection chart illustrates this by directly comparing oral fluid detection times to urine.

Drugs & Pharmaceuticals Detected

AMP	Amphetamine	50 ng/ml
MET	Methamphetamine	35 ng/ml (MDMA 100)
COC	Cocaine/Benzoylcegonine	20 ng/ml
OPI	Morphine	20 ng/ml
THC	Δ^9 -tetrahydrocannabinol	5 ng/ml
BZD	Benzodiazepines	15 ng/ml



SBL Testing Technologies Inc. (SBL) offers the leading oral fluid drug testing and breath alcohol technologies available globally. SBL also provides an unparalleled platform of comprehensive support services to ensure our client's success with these technologies. From legal and policy to training and technical, SBL's experienced team of subject matter experts will ensure success.

At SBL, we understand the sensitive nature of workplace drug & alcohol testing. Our experience and expertise will ensure that your company is ready to integrate and optimize all the advantages of oral fluid drug testing. With a corporate philosophy of unwavering commitment to service, SBL is your partner to enhance your drug & alcohol program: seamlessly and cost-effectively while successfully addressing the dynamic challenges faced by corporate drug & alcohol programs today in Canada.

In Canada, the DrugTest 5000 is for **FORENSIC USE ONLY**

403 966 5450
info@sbltt.com
www.sbltt.com

